

**OPERATE
YOUR
SYSTEM
ACCORDING
TO THIS
MANUAL
USE ONLY
THE POWER
DEVICE
SUPPLIED
WITH THIS
PRODUCT**

No user maintenance of this product is required. If the System does not function, be sure...
 (a) The TRANSFORMER or ADAPTOR is plugged into an electric outlet.
 (b) The TRANSFORMER or ADAPTOR cord is plugged into, and fully engaged in its receptacle.
 If there is a malfunction of any kind, contact your dealer for instructions on returning it for examination and repair.

MAINTENANCE INSTRUCTIONS:

1. Your Relaxor® product should never be left unattended when plugged in. Unplug from outlet when not in use, and before putting on or taking off parts.
2. Use this product only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
3. Never operate the Relaxor® product if it has a damaged cord or plug, or if it is not working properly. If it has been dropped or damaged, or dropped into water. Contact your dealer for instructions on returning it for examination or repair.
4. Close supervision is necessary when this appliance is used by, on, or near children or disabled persons.
5. Do not use outdoors.
6. Keep the cord away from heated surfaces.
7. Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
8. If you are subject to any medical condition involving the use of medical device products, and/or if you are under a health professional or physician's care, please contact them before using this product.

WARNING! To reduce risk of burns, fire, electric shock, or injury to persons.

1. Always unplug the product from the electrical outlet immediately after using and before cleaning.
2. Do not reach for a product that has fallen into water. Unplug immediately.
3. Do not use while bathing or in a shower.

DANGER! To reduce risk of electric shock.

Read all instructions before using the Relaxor® System.

IMPORTANT SAFETY INSTRUCTIONS



**Relaxor®
Ultra
Massage System**



OWNER'S MANUAL and OPERATING INSTRUCTIONS

Welcome to the world of Massage!

You are now the proud owner of another great Relaxor® product.

Your new Relaxor® massage system brings the comforting world of massage technology and places it at your fingertips. Now you can enjoy the relaxing feel of vibrating massage on your tired, knotted muscles by choosing from the dozens of pre-programed or self-programed settings.

We know you will enjoy your new Relaxor® massage system product.

**1-562-447-1780 ext. 204
FAX: 1-562-447-1798
www.inseat.com**

call or fax us at:
 11839 Smith Avenue, Santa Fe Springs, CA 90670
 Write to: Customer Service InSeat Solutions, Inc. Corporate Office
 If you need assistance for any reason, please contact us:

If there is a malfunction of any kind, contact your dealer for instructions on returning it for examination and repair.

(a) The TRANSFORMER or ADAPTOR is plugged into an electric outlet.
 (b) The TRANSFORMER or ADAPTOR cord is plugged into, and fully engaged in its receptacle.

No user maintenance of this product is required. If the System does not function, be sure...

MAINTENANCE INSTRUCTIONS:

1. Your Relaxor® product should never be left unattended when plugged in. Unplug from outlet when not in use, and before putting on or taking off parts.
2. Use this product only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
3. Never operate the Relaxor® product if it has a damaged cord or plug, or if it is not working properly. If it has been dropped or damaged, or dropped into water. Contact your dealer for instructions on returning it for examination or repair.
4. Close supervision is necessary when this appliance is used by, on, or near children or disabled persons.
5. Do not use outdoors.
6. Keep the cord away from heated surfaces.
7. Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
8. If you are subject to any medical condition involving the use of medical device products, and/or if you are under a health professional or physician's care, please contact them before using this product.

WARNING! To reduce risk of burns, fire, electric shock, or injury to persons.

1. Always unplug the product from the electrical outlet immediately after using and before cleaning.
2. Do not reach for a product that has fallen into water. Unplug immediately.
3. Do not use while bathing or in a shower.

DANGER! To reduce risk of electric shock.

Read all instructions before using the Relaxor® System.

IMPORTANT SAFETY INSTRUCTIONS

Limited Manufacturer's Warranty

Patents Pending

© Copyright 2007 InSeat Solutions, LLC. Santa Fe Springs, CA

RU85/RU85H

03/08

#111100U/#111100U-00

OPERATING INSTRUCTIONS

POWER

Activates Massage System

- Turn On: Press PWR.
- Turn Off: Press PWR again.
- 15 and 30 minute Timers:
 - Green light = 15 minute Timer.
 - Hold down Power:
 - Amber light = 30 minute Timer.
- Unit may also be turned on by pressing any ZONE 1-5 key.
- To restart after Timer shuts off:
 - Press PWR or any ZONE 1-5 key.

HEAT

Activates Heat.

- Press LO (Amber light) to turn Low Heat on.
- Press HI (Red light) to turn High Heat on.

Works with massage or independently.
30 minute Timer.

To continue after Timer shuts off, press LO or HI button again.

ZONES

(Green light indicates active Zone)

- 1 - Upper Back
- 2 - Lower Back
- 3 - Hip Region

INTENSITY

Controls INTENSITY during all modes.
(+) Increase Intensity; (-) Decrease Intensity.

SPEED

Controls pulsing frequency of massage motors in PULSE mode and controls the velocity as it cycles through in WAVE mode.

Press (+) to increase and (-) to decrease.



MODES

(Green light indicates active Mode)

ZONE -

Lets you choose any or all zones for massage.

- Press SELECT.
- Press any or all Zones (1-3) to activate.
- Press any or all Zones (1-3) to cancel.

PULSE -

Massage motors turn On and Off to simulate a pulsing action.

- Press PULSE.
- Press any or all Zones (1-3) to activate.
- Press any or all Zones (1-3) to cancel.

WAVE -

Massage motors move up and down in a wave-like motion.

- Press WAVE.
- Press any or all Zones (1-3) to activate.
- Press any or all Zones (1-3) to cancel.

TAP -

Massage rapidly turns On and Off at high Intensity to provide a tapping action.

- Press TAP.
- Press any or all Zones (1-3) to activate.
- Press any or all Zones (1-3) to cancel.

SWELL -

Massage cycle from minimum to maximum intensity to provide a swelling action

- Press SWELL.
- Press any or all Zones (1-3) to activate.
- Press any or all Zones (1-3) to cancel.

SAVE THESE INSTRUCTIONS!